


EpiGe Protocol

Adapted EpiGe Protocol to the Applied Biosystems QuantStudio 6 Flex Real-Time PCR System.

Material



- 6 rhAmp assays (IDT) for each one of the 6 cytosines of the ^{Epi}WNT-SHH panel.

Cytosine	rhAmp_ID
S1_1033	CD.GT.SPDR8897.10
S3_1292	CD.GT.GYJV9231.1
W1_2554	CD.GT.JBHB3172.1
W3_0222	CD.GT.JYVK6101.1
G1_1884	CD.GT.WFVP0601.1
G3_0126	CD.GT.PVWC1084.1

- 12 Synthetic Controls (gBlocks): Two controls (Methylated and Unmethylated) for each cytosine of the ^{Epi}WNT-SHH panel.
- rhAmp® rhAmp Genotyping Master Mix.
- rhAmp® Reporter Mix w/Reference. Use the Reporter Mix with or without reference dye, as indicated in the table below or by checking with the manufacturer of your instrument:

PCR system	Reference dye required	
	Yes	No
7900HT Fast Real-Time PCR System (Thermo Fisher Scientific)	X	
StepOne™ and StepOnePlus™ Real-Time PCR System (Thermo Fisher Scientific)	X	
Mx3005P™ and Mx4000P™ qPCR System (Agilent)	X	
7500 Real-Time PCR System (Thermo Fisher Scientific)	X	
Viia™7 Real-Time PCR System (Thermo Fisher Scientific)	X	
QuantStudio™ Flex Systems (Thermo Fisher Scientific)	X	
Biomark™ HD (Fluidigm)	X	
CFX, iQ™, and Opticon™ Real-Time PCR Detection Systems (Bio-Rad)		X
LightCycler® Real-Time PCR Systems (Roche)		X

* For instruments not listed, please check with the manufacturer.

Method

1. Use the template "EpiGe_Genotyping_Template.xlsx" to prepare the working Master mix combining the Reporter Mix with Dye and the rhAmp Master Mix in a 1:20 proportion. Take into account that all analysed samples should be replicated at least 2 times within each experiment. EpiGe-App only accepts qPCR results in which one sample has been analysed at the time and have at least 2 NTC replicates for each assay.
2. To a new tube or vial, add the following:

	Volume
rhAMP master Mix	100 ul
rhAMP Reporter Mix with Dye	5ul

3. Mix and vortex.
4. Prepare the six EpiGe SNP Genotyping Assay Reaction Mixes (Combined Master Mix, Reporter Mix, and EpiGe assay) to a final volume of 5µl for 6 samples (2 NTC, 2 Samples and 2 additional reactions to account for pipetting errors).

	Samples	W1_2554	W3_0222	S1_1033	S3_1292	G1_1884	G3_0126
	5ul	6	6	6	6	6	6
	Final	Samples	# of Samples	# of Samples	# of Samples	# of Samples	# of Samples
Master Mix + Reporter	2.65	15.9	15.9	15.9	15.9	15.9	15.9
SNP assay (20X)	0.25	1.5	1.5	1.5	1.5	1.5	1.5
Water	1.1	6.6	6.6	6.6	6.6	6.6	6.6
Mix Volume	4	24	24	24	24	24	24

Sample	1	1
Final Volume	5	5

5. Vortex and briefly centrifuge before use.
6. Add 4µl of EpiGe Master Mix to each well of the qPCR plate or strip and 1µl of the bisulfite converted sample or DNA synthetic control (gBlock).

	1	2	3	4	5	6	7	8	9	10	11	12
A				W1_2554 NTC	W3_0222 NTC	S1_1033 NTC	S3_1292 NTC	G1_1884 NTC	G3_0126 NTC			
B				W1_2554 NTC	W3_0222 NTC	S1_1033 NTC	S3_1292 NTC	G1_1884 NTC	G3_0126 NTC			
C				W1_2554 Sample	W3_0222 Sample	S1_1033 Sample	S3_1292 Sample	G1_1884 Sample	G3_0126 Sample			
D				W1_2554 Sample	W3_0222 Sample	S1_1033 Sample	S3_1292 Sample	G1_1884 Sample	G3_0126 Sample			
E				W1_2554 gBlock_W1_U	W3_0222 gBlock_W3_U	S1_1033 gBlock_S1_U	S3_1292 gBlock_S3_U	G1_1884 gBlock_G1_U	G3_0126 gBlock_G3_U			
F				W1_2554 gBlock_W1_U	W3_0222 gBlock_W3_U	S1_1033 gBlock_S1_U	S3_1292 gBlock_S3_U	G1_1884 gBlock_G1_U	G3_0126 gBlock_G3_U			
G				W1_2554 gBlock_W1_M	W3_0222 gBlock_W3_M	S1_1033 gBlock_S1_M	S3_1292 gBlock_S3_M	G1_1884 gBlock_G1_M	G3_0126 gBlock_G3_M			
H				W1_2554 gBlock_W1_M	W3_0222 gBlock_W3_M	S1_1033 gBlock_S1_M	S3_1292 gBlock_S3_M	G1_1884 gBlock_G1_M	G3_0126 gBlock_G3_M			