A) Combined MasterMix and Reporter Mix

Final Volume (Multiple rhAMP master Mix 1000ul rh Reporter Mix with 50ul Vol **100** ul 5 ul

B) Samples and Primers

	Reporter	Quencher	Task
Allele1 (Meth)	FAM	NFQ-MGB	1/1
Allele2 (Unmeth)	VIC	NFQ-MGB	2/2

Primers

Primer Name	ID				
W1_2554	CD.GT.JBHB3172.1				
W3_0222	CD.GT.JYVK6101.1				
S1_1033	CD.GT.SPDR8897.10				
S3_1292	CD.GT.GYJV9231.1				
G1_1884	CD.GT.WFVP0601.1				
G3_0126	CD.GT.PVWC1084.1				

C) Working Mix

	Samples	W1_2554	W3_0222	S1_1033	S3_1292	G1_1884	G3_0126
	5ul	6	6	6	6	6	6
	Final	Samples	# of Samples				
Master Mix + Reporter	2.65	15.9	15.9	15.9	15.9	15.9	15.9
SNP assay (20X)	0.25	1.5	1.5	1.5	1.5	1.5	1.5
Water	1.1	6.6	6.6	6.6	6.6	6.6	6.6
Mix Volume	4	24	24	24	24	24	24

Sample	1	1
Final Volume	5	5

D) qPCR Template

	1	2	3	4	5	6	7	8	9	10	11	12
^				W1_2554	W3_0222	S1_1033	S3_1292	G1_1884	G3_0126			
Α				NTC	NTC	NTC	NTC	NTC	NTC			
В				W1_2554	W3_0222	S1_1033	S3_1292	G1_1884	G3_0126			
				NTC	NTC	NTC	NTC	NTC	NTC			
•				W1_2554	W3_0222	S1_1033	S3_1292	G1_1884	G3_0126			
·				Sample	Sample	Sample	Sample	Sample	Sample			
D				W1_2554	W3_0222	S1_1033	S3_1292	G1_1884	G3_0126			
				Sample	Sample	Sample	Sample	Sample	Sample			
E				W1_2554	W3_0222	S1_1033	S3_1292	G1_1884	G3_0126			
-				gBlock_W1_U	gBlock_W3_U	gBlock_S1_U	gBlock_S3_U	gBlock_G1_U	gBlock_G3_U			
				W1_2554	W3_0222	S1_1033	S3_1292	G1_1884	G3_0126			
•				gBlock_W1_U	gBlock_W3_U	gBlock_S1_U	gBlock_S3_U	gBlock_G1_U	gBlock_G3_U			
G				W1_2554	W3_0222	S1_1033	S3_1292	G1_1884	G3_0126			
•				gBlock_W1_M	gBlock_W3_M	gBlock_S1_M	gBlock_S3_M	gBlock_G1_M	gBlock_G3_M			
н				W1_2554	W3_0222	S1_1033	S3_1292	G1_1884	G3_0126			
••				gBlock_W1_M	gBlock_W3_M	gBlock_S1_M	gBlock_S3_M	gBlock_G1_M	gBlock_G3_M			

E) Thermal Cycling

Step	Temperature (°C)	Time (min:sec)	Cycles
Enzyme Activation	95	10:00	1
Denautre	95	0:10	
Anneal	60	0:30	40
Extend	68	0:45	
PCR- read	60	0:30	1